# **ZEE nit 700**

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

The image shows a software window titled "Integration parameters" with several tabs: "Optical parameters", "Energy / Gain", "Integration parameters" (selected), "Scan parameters", and "Scan". The window contains four main sections:

- Integration mode:** Radio buttons for "Mean", "Repeated mean" (selected), "Peak area", and "Peak height".
- Measuring times:** Input fields for "Integration:" (3.0 s), "Autozero:" (3.0 s), and "Delay:" (5.0 s).
- Autozero mode:** Radio buttons for "Simple mode", "Wash solution" (selected), and "AZ solution". A checkbox "Disable AZ drift check" is at the bottom.
- Peak smoothing:** Spinners for "Average:" (20 Points) and "Golay-Savitzky:" (25 Points).

At the bottom right, there are two buttons: "Opt.meas.times" and "Reset".

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
N2O: OK Gas/Ox: l/h

Test air Test N2O Test fuel Test end

Ignite flame Air -> N2O Extinguish flame

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

The screenshot shows the 'Autosampler' configuration window with several tabs: 'Autosample mode', 'Positions', 'Techn. parameters', 'Dilution/IB', 'Function tests', and 'Error check'. The 'Autosample mode' tab is active. It contains three main sections: 'Autosampler type' with radio buttons for 'No autosampler', 'AS51', 'AS52', 'AS52 S', 'MPE50', and 'MPE60'; 'Tray type' with radio buttons for '89 Positions', '53 Positions', and '50 Positions'; and 'Working mode' with radio buttons for 'Step [1 cycle]' and 'Continuous'. The 'Washing' section has radio buttons for 'Off', 'Between samples', 'Between measurements', and 'Between components', along with input fields for 'Wash cycles' and 'Wash time [s]'. The 'Injection switch' section has radio buttons for 'Off' and 'Activate SFS 5', and input fields for 'Load time' and 'Injection time'. The 'Sample volume' section has input fields for 'Injection volume' and 'Dip time'. A 'Reset' button is located at the bottom right of the main configuration area, and an 'OK' button is at the bottom right of the window.

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The screenshot shows the 'Calibration' window with tabs: 'Calibration mode', 'Conditions', 'Statistics', 'Conc. input', 'Table', 'Curve parameters', and 'Recalibration'. The 'Calibration mode' tab is active. It contains a section titled 'Computation / calibration method:' with radio buttons for 'Absorbance mode (no calibration)', 'Standard calibration', 'Bracketing calibration', 'Method of addition', and 'Method of Additions calibration'. A 'Reset' button is located at the bottom right of the main configuration area, and an 'OK' button is at the bottom right of the window.

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The screenshot shows the 'Conditions' tab of the Calibration software. The interface includes several sections for setting up the calibration process:

- Standards:**
  - Making standards:** Radio buttons for 'Using premixed standards' (selected), 'By different volumes', and 'Mixed by autosampler'.
  - No. of standards:** A numeric input field set to 5.
  - Sample meas. without add:** A numeric input field set to 3.
- Blank correction:** Radio buttons for 'Absorbance corrected' (selected) and 'Concentration corrected'.
- Other options:** Checkboxes for 'Continue after curve calculation' and 'AZ (AZ/AM) between standards'.
- Getting values by:** Radio buttons for 'Measurement' (selected) and 'Manual input'.
- Mix conditions:**
  - Filling level mix cup:** A numeric input field set to 5 ml.
  - [sample part]:** A numeric input field set to 2.5 ml.
  - Sample volume:** A numeric input field set to 0 ml.
  - Max. add. volume:** A numeric input field set to 10 ml.
  - Total modifier vol.:** A numeric input field set to 0 ml.

Buttons for 'Reset', 'OK', and 'Recalibration' are visible at the bottom of the window.

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The screenshot shows the 'Statistics' tab of the Calibration software. The interface includes several sections for setting up statistical analysis:

- Statistics:** Radio buttons for 'Off', 'Sigma statistics' (selected), and 'Median Statistics'.
- Grubbs outlier test:** Radio buttons for 'Off' (selected), 'On (Flagging mean by "[I]")', and 'On [Without flagging]'.
- Measurement cycles:**
  - Runs/Std.:** A numeric input field set to 3.
  - Blind runs/Std.:** A numeric input field set to 1.
- Statistical output:** Checkboxes for 'Standard deviation [SD]' and 'Relat. standard deviation [RSD]' (both selected).

Buttons for 'Reset' and 'OK' are visible at the bottom of the window.

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☐ Solid

Concentrations:

|                                            | Conc.  |
|--------------------------------------------|--------|
| Stock solution 1:                          | 1.0000 |
| <input type="checkbox"/> Stock solution 2: | 1.0000 |
| <input type="checkbox"/> Stock solution 3: | 1.0000 |
| <input type="checkbox"/> Stock solution 4: | 1.0000 |

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
Clear Table Set recal Check Run sample

Delete line Parameters Fit curve -> Blindw.

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
The "Standard parameters" dialog box contains the following fields and options:

- Standard:** A dropdown menu set to "Cal-Std1".
- Sampler position:** A text box containing "0" and a small unit icon.
- Making standards:** A section containing:
  - Mixing ratio:** A text box with "0" and a percentage sign.
  - from Stock std.:** A group box containing four radio buttons:
    - ☒ #1 (1,000 mg/L)
    - ☐ #2 (1,000 mg/L)
    - ☐ #3 (1,000 mg/L)
    - ☐ #4 (1,000 mg/L)
- Concentration:** A text box with "0.1000" and a unit "mg/L".
- Absorbance:** A text box with "Mean" below it.
- Buttons:** "Cancel" and "OK" at the bottom right.

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The "Calibration" dialog box has several tabs: "Calibration mode", "Conditions", "Statistics", "Conc. input", "Table", "Curve parameters" (selected), and "Recalibration". The "Curve parameters" tab contains the following sections:

- Curve calculation:**
  - Curve form:** Radio buttons for "Linear", "Nonlinear", and "Automatically" (selected).
  - Intercept:** Radio buttons for "Zero" and "Calculate" (selected).
- Weights of calibration points:** Radio buttons for "No weights" (selected), "1 / Conc", "1 / SD", and "1 / (SD \* Conc)".
- Curve examination:** Radio buttons for "No outlier test" (selected), "Flag and eliminate outlier", and "Repeat measurement of the outlier".
- Buttons:** "Reset" and "OK" at the bottom right.

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





**Samples**

Sample table | Conc. output | ID samples | Statistics | Mult.Param.

| No | * | Solution |     | ID/Wt-Data |        |          | Actions |       |      |
|----|---|----------|-----|------------|--------|----------|---------|-------|------|
|    |   | Type     | Pos | Name       | Pre-DF | Vol.[mL] | Wt[g]   | Cal.  | Tube |
| 1  | * | Sample   | ### |            | 1.000  |          |         | Rec/C |      |
| 2  | * | Sample   | ### |            | 1.000  |          |         |       |      |
| 3  | * | Sample   | ### |            | 1.000  |          |         |       |      |
| 4  | * | Sample   | ### |            | 1.000  |          |         |       |      |
| 5  | * | Sample   | ### |            | 1.000  |          |         |       |      |
| 6  |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 7  |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 8  |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 9  |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 10 |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 11 |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 12 |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 13 |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 14 |   | Sample   | ### |            | 1.000  |          |         |       |      |
| 15 |   | Sample   | ### |            | 1.000  |          |         |       |      |

Buttons: Init.Table, Change, Load / Save, Single values, Final actions (Lamp off, Flame off), Delete results, Delete, Working area, -> Cal.tab., Print table, Insert, Run sample, Show sequence, Recalculate, Export [CSV], Clipboard, OK

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Sample table → ID    Change line    Insert line    Save    → Sample table

Delete table    Delete    Append line    Load    Synopsis

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**Samples**

Sample table | Conc. output | ID samples | **Statistics** | Mult. Param. |

Statistics:

☐ Off  
☒ Sigma statistics  
☐ Median Statistics

Cycles:

Measurement cycles:   
Runs/QC sample:   
Blind cycles:

Statistical output:

☒ Standard deviation (SD)  
☒ Relat. standard deviation (RSD)  
☐ Relative homogeneous factor (RHF)

Grubbs outlier test:

☒ Off  
☐ On (Flagging mean by "[I]")  
☐ On (Without flagging)

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**Samples**

Sample table | Conc. output | ID samples | Statistics | Mult.Param.

| No | * | Solution |     | ID/Wt.Data |        |          |       | Actions |      |    |
|----|---|----------|-----|------------|--------|----------|-------|---------|------|----|
|    |   | Type     | Pos | Name       | Pre-DF | Vol.[mL] | Wt[g] | Cal.    | Tube | AZ |
| 1  | * | Sample   | ### |            | 1.000  |          |       | Rec/C   |      |    |
| 2  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 3  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 4  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 5  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 6  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 7  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 8  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 9  | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 10 | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 11 | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 12 | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 13 | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 14 | * | Sample   | ### |            | 1.000  |          |       |         |      |    |
| 15 | * | Sample   | ### |            | 1.000  |          |       |         |      |    |

Buttons: Init.Table, Change, Load / Save, Single values, Final actions, Delete results, Delete, Working area, -> Cal.tab., Lamp off, Flame off, Print table, Insert, Run sample, Show sequence, Recalculate, Export [CSV], Clipboard, OK

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The screenshot shows the 'Autosampler' window with the 'Specials' tab active. The 'Modifiers' section has three rows for #1, #2, and #3, each with a 'Name' field and a 'Vol. [μl]' spinner set to 10. The 'Enrichment in tube' section has a 'Mode' group with radio buttons for 'Off' (selected), 'Permanent [Samples only]', 'If conc. too low', and 'Permanent [incl. calib.]'. Below this is an 'Enrichment cycles' spinner set to 2 and a checkbox for 'Repeat before enrichment'. The 'Thermal pretreatment' section has an 'Activation' group with 'Off' (selected) and 'On'. The 'Warm up components' section has radio buttons for 'Sample' (selected), 'Modifier 1', 'Sample + Modifier 1', 'Modifier 1 + Sample', and 'Modifier 1 + Modifier 2'. A 'Warm up delay' spinner is set to 0. A 'Reset' button is at the bottom right of the main panel, and an 'OK' button is at the very bottom.

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The screenshot shows the 'Calibration' window with the 'Calibration mode' tab active. The 'Computation / calibration method' section has radio buttons for 'Absorbance mode (no calibration)', 'Standard calibration' (selected), 'Bracketing calibration', 'Method of addition', and 'Method of Additions calibration'. A 'Reset' button is located at the bottom right of the main panel, and an 'OK' button is at the very bottom.

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The screenshot shows the 'Calibration' window with the 'Conditions' tab active. The interface includes the following sections:

- Standards:**
  - Making standards:** Radio buttons for 'Using premixed standards', 'By different volumes', and 'By dilution' (selected).
  - No. of standards:** A numeric input field set to 5.
  - Sample meas. without add:** A numeric input field set to 3.
- Blank correction:**
  - Radio buttons for 'Absorbance corrected' (selected) and 'Concentration corrected'.
  - Checkboxes for 'Continue after curve calculation' and 'AZ [AZ/AM] between standards'.
- Getting values by:** Radio buttons for 'Measurement' (selected) and 'Manual input'.
- Mix conditions:**
  - Filling level mix cup:** A numeric input field set to 5 ml.
  - [sample part]:** A numeric input field set to 2.5 ml.
  - Sample volume:** A numeric input field set to 20 ml.
  - Max. add. volume:** A numeric input field set to 10 ml.
  - Total modifier vol.:** A numeric input field set to 0 ml.

Buttons for 'Reset', 'OK', and 'Cancel' are located at the bottom right.

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The screenshot shows the 'Calibration' window with the 'Statistics' tab active. The interface includes the following sections:

- Statistics:** Radio buttons for 'Off', 'Sigma statistics' (selected), and 'Median Statistics'.
- Measurement cycles:**
  - Runs/Std.:** A numeric input field set to 3.
  - Blind runs/Std.:** A numeric input field set to 1.
- Grubbs outlier test:** Radio buttons for 'Off' (selected), 'On (Flagging mean by "[I]")', and 'On [Without flagging]'.
- Statistical output:** Checkboxes for 'Standard deviation [SD]' and 'Relat. standard deviation [RSD]' (both selected).

Buttons for 'Reset', 'OK', and 'Cancel' are located at the bottom right.

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☒ Liquid  
☐ Solid

Concentrations:

|                                            | Conc.   |
|--------------------------------------------|---------|
| Stock solution 1:                          | 2.50000 |
| <input type="checkbox"/> Stock solution 2: | 1.0000  |
| <input type="checkbox"/> Stock solution 3: | 1.0000  |
| <input type="checkbox"/> Stock solution 4: | 1.0000  |

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☒ No weights  
☐ 1 / Conc  
☐ 1 / SD  
☐ 1 / (SD \* Conc)

**Curve examination:**

☒ No outlier test  
☐ Flag and eliminate outlier  
☐ Repeat measurement of the outlier

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Clear Table | Set recal | Check | Run sample

Delete line | Parameters | Fit curve | → Blindw.

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


Sample table → ID | Change line | Insert line | Save | → Sample table

Delete table | Delete | Append line | Load | Synopsis

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